

Product Datasheet

Human Immunodeficiency Virus 2 (gp36) Antigen, Recombinant V9003

Product Name Human Immunodeficiency Virus 2 (gp36) Antigen, Recombinant

Catalog Number V9003

Species/Host E. coli

Tested Applications LF, ELISA, WB

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note His tag

Purity ~95% by SDS-PAGE

Buffer 20 mM Carbonate

Specificity Ectodomain, gp36

Molecular Weight 35.8 kDa