

Product Datasheet

Treponema pallidum (TnpA) Antigen, Recombinant V9004

Product Name Treponema pallidum (TnpA) Antigen, Recombinant

Catalog Number V9004

Species/Host E. coli

Tested Applications LF, ELISA, WB

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note His tag

Purity ~95% by SDS-PAGE

Buffer 20 mM Na-carbonate, pH 10.0

Specificity TnpA

Molecular Weight 16.4 kDa