

Product Datasheet

Treponema pallidum (TpN47) Antigen, Recombinant V9006

Product Name Treponema pallidum (TpN47) Antigen, Recombinant

Catalog Number V9006

Species/Host E. coli

Tested Applications LF, ELISA, WB

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note His tag

Purity ~95% by SDS-PAGE

Buffer 20mM Carbonate

Specificity TpN47

Molecular Weight 47.7 kDa