

Product Datasheet

Hepatitis C virus (core) Antigen, Recombinant V9007

Product Name Hepatitis C virus (core) Antigen, Recombinant

Catalog Number V9007

Species/Host E. coli

Tested Applications LF, ELISA, WB

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage $\leq -20^{\circ}\text{C}$

Note His tag

Purity ~95% by SDS-PAGE

Buffer Phosphate buffered saline pH 7.4

Specificity Core protein

Molecular Weight 48.7 kDa