

## Product Datasheet

### Hepatitis C virus (core) Antigen, Recombinant V9007

**Product Name** Hepatitis C virus (core) Antigen, Recombinant

**Catalog Number** V9007

**Species/Host** E. coli

**Tested Applications** LF, ELISA, WB

**Brand** ViroStat

**Form/Appearance** Purified recombinant antigen

**Storage** ≤ -20°C

**Note** His tag

**Purity** ~95% by SDS-PAGE

**Buffer** Phosphate buffered saline pH 7.4

**Specificity** Core protein

**Molecular Weight** 48.7 kDa