

Product Datasheet

Anti-Human CD21 Monoclonal Antibody, Unconjugated Clone B-E5 854.163.020

Product Name Anti-Human CD21 Monoclonal Antibody, Unconjugated Clone B-E5

Catalog Number 854.163.020

Description Mouse monoclonal

Species/Host Mouse

Conjugation Purified

Tested Applications IHC, FC

Immunogen Daudi cell line

Brand Diaclone

Form/Appearance Phosphate-buffered saline with 1% BSA and 0.09% sodium azide

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Note Use 10 µl to label 10⁶ cells or 100 µl of whole blood

Isotype Mouse IgG2a Kappa light chain

Clonality Monoclonal

Specificity Recognises the Epstein-Barr virus receptor (EBV R), CR2, C3d R, a 110 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells