

Product Datasheet

Herpes Simplex Virus 1 (gG) Antigen, Recombinant V9010

Product Name	Herpes Simplex Virus 1 (gG) Antigen, Recombinant
Catalog Number	V9010
Species/Host	E. coli
Tested Applications	ELISA, WB
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	C-terminus 6-his; Both HSV-1 and HSV-2 positive human sera are reactive with this antigen. Reactive with V0198, V0199, V0194, V0190, & V0192.
Purity	~95% by SDS-PAGE
Buffer	PBS
Specificity	Glycoprotein gD, a.a. 226 thru a.a. 394