

## Product Datasheet

### Herpes Simplex Virus 2 (gD) Antigen, Recombinant V9011

**Product Name** Herpes Simplex Virus 2 (gD) Antigen, Recombinant

**Catalog Number** V9011

**Species/Host** E. coli

**Tested Applications** ELISA, WB

**Brand** ViroStat

**Form/Appearance** Purified recombinant antigen

**Storage** ≤ -20°C

**Note** No fusion partner present

**Purity** ~95% by SDS-PAGE

**Buffer** PBS

**Specificity** Glycoprotein gD, a.a. 31 thru a.a. 335

**Molecular Weight** 62.8 kDa