

Product Datasheet

Herpes Simplex Virus 2 (gD) Antigen, Recombinant V9011

Product Name Herpes Simplex Virus 2 (gD) Antigen, Recombinant

Catalog Number V9011

Species/Host E. coli

Tested Applications ELISA, WB

Brand ViroStat

Form/Appearance Purified recombinant antigen

Storage ≤ -20°C

Note No fusion partner present

Purity ~95% by SDS-PAGE

Buffer PBS

Specificity Glycoprotein gD, a.a. 31 thru a.a. 335

Molecular Weight 62.8 kDa