

Product Datasheet

Cytomegalovirus (pp65) Antigen, Recombinant V9012

Product Name	Cytomegalovirus (pp65) Antigen, Recombinant
Catalog Number	V9012
Species/Host	E. coli
Tested Applications	ELISA
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Note	GST
Purity	~95% by SDS-PAGE
Buffer	25 mM Tris-HCl, pH 7.2 + 1 mM EDTA + 50% glycerol
Specificity	pp65 (UL83), 561 a.a.