

Product Datasheet

Cytomegalovirus (pp65) Antigen, Recombinant V9012

| | |
|----------------------------|---|
| Product Name | Cytomegalovirus (pp65) Antigen, Recombinant |
| Catalog Number | V9012 |
| Species/Host | E. coli |
| Tested Applications | ELISA |
| Brand | ViroStat |
| Form/Appearance | Purified recombinant antigen |
| Storage | ≤ -20°C |
| Note | GST |
| Purity | ~95% by SDS-PAGE |
| Buffer | 25 mM Tris-HCl, pH 7.2 + 1 mM EDTA + 50% glycerol |
| Specificity | pp65 (UL83), 561 a.a. |