

Product Datasheet

Cytomegalovirus (pp150) Antigen, Recombinant V9013

| | |
|----------------------------|--|
| Product Name | Cytomegalovirus (pp150) Antigen, Recombinant |
| Catalog Number | V9013 |
| Tested Applications | ELISA, WB |
| Brand | ViroStat |
| Form/Appearance | Purified recombinant antigen |
| Storage | ≤ -20°C |
| Purity | ~95% by SDS-PAGE |
| Buffer | 25 mM Tris, pH 8.0 + 1.5 M Urea + 50% glycerol + 0.2% Triton |
| Specificity | pp150 (UL32), a.a.426 to a.a.606 |