

Product Datasheet

Cytomegalovirus (pp150) Antigen, Recombinant V9013

Product Name	Cytomegalovirus (pp150) Antigen, Recombinant
Catalog Number	V9013
Tested Applications	ELISA, WB
Brand	ViroStat
Form/Appearance	Purified recombinant antigen
Storage	≤ -20°C
Purity	~95% by SDS-PAGE
Buffer	25 mM Tris, pH 8.0 + 1.5 M Urea + 50% glycerol + 0.2% Triton
Specificity	pp150 (UL32), a.a.426 to a.a.606