

Product Datasheet

Oropouche Virus (OROV-NP) Antibody HM1405

Product Name	Oropouche Virus (OROV-NP) Antibody
Catalog Number	HM1405
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant Oropouche virus
Alternative Names	OPV, OPOV
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% NaCl
Specificity	Oropouche virus nucleoprotein (NP)
Matched Pairs	HM1407, HM1406 or HM1405 (detection)
Preservatives	0.05% Sodium Azide