

MDMA-BSA Antigen

Catalog number	170056
Description	Synthetic MDMA, BSA conjugated
Amino acid sequence	-
Product host	-
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	60 months from manufacturing at -20 °C. Short term storage (≤ 2 weeks): 2–8 °C. Avoid repeated freezing and thawing.
Analyte description	MDMA is a synthetic drug also known as Ecstasy. Usage may alter mood, heart rate, blood pressure, activity and energy levels, and lead to visual and auditory hallucinations.
Product concentration	Lot-specific
Matched pairs	100906 (detection)
Miscellaneous	This product was previously available with a Catalog code P82-99-51A. It has been produced using the same conjugation protocol.

**Legal disclaimer**