

## Amphetamine-HRP Antigen

---

<b>Catalog number</b>	170100
<b>Description</b>	Synthetic Amphetamine, HRP conjugated
<b>Amino acid sequence</b>	-
<b>Product host</b>	-
<b>Product formulation</b>	Liquid; Antigen-HRP conjugate
<b>Product buffer solution</b>	1:1 0.1 M Tris base / FBS, pH 7.0, 0.0016 % Gentamicin sulfate, 0.005 % thimerosal
<b>Reconstitution</b>	-
<b>Shelf life and storage</b>	Unspecified, storage at 2–8 °C. Do not freeze.
<b>Analyte description</b>	Amphetamine is a central nervous system stimulant drug used to treat ADHD
<b>Product concentration</b>	-
<b>Purity</b>	-
<b>Matched pair</b>	C03-99-72M-P (capture)
<b>Miscellaneous</b>	This product was previously available with Catalog code P91-99-71H. It has been produced using the same conjugation protocol. Linked via para position.



**Legal disclaimer**