

## Testosterone-3-HRP Antigen

---

<b>Catalog number</b>	170068
<b>Description</b>	Synthetic Testosterone-3, HRP conjugated
<b>Amino acid sequence</b>	-
<b>Product host</b>	-
<b>Product formulation</b>	Liquid; Antigen-HRP conjugate
<b>Product buffer solution</b>	1:1 0.1 M Tris base / FBS, pH 7.0, 0.0016% Gentamicin sulfate, 0.005% thimerosal
<b>Reconstitution</b>	-
<b>Shelf life and storage</b>	Shelf life unspecified, storage at 2–8 °C. Do not freeze.
<b>Analyte description</b>	Testosterone is a hormone produced by reproductive organs. Higher amounts are found in males and responsible for the development of male sex characteristics. Testosterone levels can be measured in blood and may affect the health of men and women.
<b>Product concentration</b>	-
<b>Purity</b>	-
<b>Antibodies tested</b>	-
<b>Miscellaneous</b>	This product was previously available with catalog code P91-92-71H. It has been produced using the same conjugation protocol.

**Legal disclaimer**