

Amphetamine-BSA Antigen

Catalog number	170069
Description	Synthetic Amphetamine, BSA conjugated
Amino acid sequence	-
Product host	-
Product formulation	Liquid; Antigen-BSA conjugate
Product buffer solution	Distilled water
Reconstitution	-
Shelf life and storage	60 months from manufacturing at -20 °C. Short term storage (≤ 2 weeks): 2–8 °C. Avoid repeated freezing and thawing.
Analyte description	Amphetamine is a central nervous system stimulant drug used to treat ADHD and narcolepsy.
Product concentration	Lot-specific
Matched pair	C03-99-72M-P (detection)
Miscellaneous	This product was previously available with a Catalog code P81-99-71A. It has been produced using the same conjugation protocol.



Legal disclaimer